

APESIN rapid plus



Liquid fast-acting Surface disinfectant cleaner

- Broad activity
- Versatile
- 2-in-1

Product profile

- APESIN rapid plus is a liquid fast-acting surface disinfectant cleaner for entire institutional, medical and food-processing areas.
- With a broad activity against bacteria, yeast, fungi and viruses.
- For simultaneous disinfecting and cleaning of washable surfaces and items.
- APESIN rapid plus is free of fragrances.
- The short contact times enable efficient work processes.

Area of application

- APESIN rapid plus should be applied with wet wiping method (mechanical action).
- Always pay attention to complete wetting of the surface and do not dry surface before completing required contact time. Follow the specified contact times for the respective product concentrations (see page 2).
- Surfaces with food contact need to be rinsed with clear water after use and contact time.
- APESIN rapid plus is suitable for surfaces and items made of e.g., PA/POM, acrylic glass, VA steel, etc. Before applying, a material compatibility examination in an inconspicuous place is recommended.
- Not applicable on metals such as brass, copper or aluminum.

Ingredients:

Ingredients per 100g: 2,475g DIDECYLDIMONIUM CHLORIDE, <5% Non-ionic surfactants

Sustainable development and production:




Product safety, storage and environmental instructions

Safety: For professional use only. Test material compatibility in an inconspicuous place before applying. For further information, see Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

Storage: Store at room temperature in the original container. Keep from freezing.

Environmental: Only give completely emptied packaging to special waste disposer.

Application and dosing



Maintenance cleaning: Mop the floor with a clean mop.



Surface cleaning: Wipe the surface wet with a clean cloth.



Sales units:

Order No.: 716035 2 x 5 L
Order No.: 716033 4 x 2 L

Your trusted partner on site:

pH-Value:



APESIN rapid plus



Liquid fast-acting Surface disinfectant cleaner

Surface Disinfection:		contact time				
Food area / industrial & institutional area		1 min.	5 min.	15 min.	30 min.	60 min.
bactericidal activity						
EN 1276	high organic load	1%	1%			1%
EN 13697	high organic load	3%	1,5%			0,75%
yeastcidal activity						
EN 1650	high organic load	0,75%	0,5%	0,5%		0,5%
EN 13697	high organic load	2%	1%	0,1%		0,25%
fungicidal activity						
EN 1650	high organic load			4%		2%
EN 13697	high organic load			3%		2%
Surface Disinfection:		contact time				
Medical area		1 min.	5 min.	15 min.	30 min.	60 min.
bactericidal activity						
EN 13727	high organic load	0,5%	0,5%	0,125%		
EN 16615	high organic load	2,5%	1%	0,25%		
yeastcidal activity						
EN 13624	high organic load	2,5%	0,5%	0,25%		
EN 16615	high organic load	2,5%	1%	0,25%		
fungicidal activity						
EN 13624	high organic load			5%		3%
		contact time				
		1 min.	5 min.	15 min.	30 min.	60 min.
virucidal activity						
EN 14476 (full spectrum)	low organic load			5%	2,5%	
EEN 14476 limited spectrum	high organic load				2,5%	
virucidal activity						
EN 14476 virucidal activity against enveloped viruses (incl. HIV, Coronavirus, Hepatitis B and C virus and human and animal influenzaviruses)	high organic load		2%			
EN 14476 individual test organisms						
EN 14476 Vacciniavirus	low organic load					
	high organic load		2%			
EN 14476 Adenovirus	low organic load			2%		
	high organic load				2,5%	
EN 14476 Norovirus	low organic load			2%		
	high organic load			2,5%		
EN 14476 Poliovirus	low organic load			5%	2,5%	
VAH-Certificate		contact time				
surface disinfection in medical areas bactericidal and yeasticidal activity		1 min.	5 min.	15 min.	30 min.	60 min.
	with mechanical action dirty conditions	2,5%	1%	0,5%	0,25%	

