


# APESIN rapid



## Liquid quick acting disinfecting surface cleaner

- Economic
- Very efficient
- Dermatologically tested

### Product Profile

- APESIN rapid is a surface disinfecting cleaner with a broad activity at a very low application concentration thanks to special combination of active ingredients.
- Preventive against skin/foot pathogens.
- APESIN rapid has an excellent cleaning action together with high material compatibility. It has a good sensory acceptance among staff using the product and patients.
- It can be used in remote metering equipment (BAM-tested = Federal Institute of Materials Testing and Research). When using the 2L dosing bottle, the integral dosing device makes metering simple and reliable. It is highly cost-effective and extremely economical because of the low concentration used, together with the fast.

### Application & Dosing



Dose in relation to application and degree of soil. Follow below instructions. See exposure time stated on the next page.



**Manual floor cleaning:** Wipe the floor with a clean mop.



**Surface cleaning:** Wipe the surface wet with a clean cloth.



### Area of application

- APESIN rapid can be used on all water resistant materials throughout the hospital sectors, old people's homes and in food-processing areas.
- Suitable for surfaces and items such as PVC / PP / PE, acrylic glass, rubber, synthetic resins, linoleum, brass, stainless steel, chrome-plated fittings, etc.
- It is also ideally suited for areas that need to be particularly well protected against infections or from where infections are likely to originate.
- Always ensure complete wetting of the surfaces.

**Ingredients:** Ingredients per 100g: 11g DIDECYLDIMONIUM CHLORIDE, 10g LAURYLAMINE DIPROPYLENEDIAMINE, 5 - <15% Cationic surfactants, Non-ionic surfactants, <5% Phosphonates, Perfumes, LIMONENE, LINALOOL, alpha-ISOMETHYL IONONE, GERANIOL, BUTYLPHENYL METHYLPROPIONAL

### Sustainable development and production:



### Product safety, storage and environmental instructions

**Safety:** For professional use only. For further information, see Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

**Storage:** Store at room temperature in the original container.

**Environmental:** Only give completely emptied container to special waste disposer.

#### Sales units:


Order-No.: 713179	6 x 2 L
Order-No.: 713178	2 x 5 L
Order-No.: 713181	4 x 2 L

#### Your trusted partner on site:

pH-value:



# APESIN rapid



## Liquid quick acting disinfecting surface cleaner


<b>Surface Disinfection:</b>		<b>contact time</b>					
Food area / industrial & institutional area		5 min.	15 min.	30 min.	60 min.	120 min.	240 min.
<b>bactericidal activity</b>							
EN 1276	low organic load	0,1%					
	high organic load	0,25%					
EN 13697	low organic load	0,25%					
	high organic load	0,50%					
<b>yeasticidal activity</b>							
EN 1650	low organic load	0,1%					
	high organic load	0,25%					
EN 13697	low organic load	0,1%					
	high organic load	0,25%					
<b>Surface Disinfection:</b>		<b>contact time</b>					
Medical Area		5 min.	15 min.	30 min.	60 min.	120 min.	240 min.
<b>bactericidal + yeasticidal activity</b>							
EN 16615	high organic load		0,5 %		0,25 %		
<b>activity against specific organisms: dermatophytes (skin fungi)</b>							
EN 13624 (Trichophyton interdigitale)	high organic load			0,25 %			
<b>additional activities and tests</b>							
tuberculocidal (Mycobacter terrae)**				3 %	1 %		
Limited virucidal against enveloped viruses* (incl. HBV, HCV, HIV, human and animal influenza viruses)			0,5%				
Adenovirus*				0,5%	0,25%		
Vaccinia Virus*			0,5%				
BVDV*		0,25%					
Rotavirus*		0,25%					
Norovirus*						1%	
Polyomavirus SV 40*							0,5%
<b>VAH -Certificate</b>		<b>contact time</b>					
Surface disinfection in medical areas		5 min.	15 min.	30 min.	60 min.	120 min.	240 min.
bactericidal and yeasticidal activity	with mechanical action dirty conditions		0,5 %		0,25 %		

\*tested according DVV-method

\*\* DGHM-method



# APESIN rapid



## Liquid quick acting disinfecting surface cleaner

<b>Surface Disinfection:</b>		<b>contact time</b>					
Food area / industrial & institutional area		5 min.	15 min.	30 min.	60 min.	120 min.	240 min.
bakterizide + levurozide Wirkung							
EN 1276 + EN 1650 + EN 13697	low organic load			0,25%			
	high organic load			0,5%			
<b>Surface Disinfection:</b>		<b>contact time</b>					
Medical area		5 min.	15 min.	30 min.	60 min.	120 min.	240 min.
bactericidal + yeasticidal activity (with mechanical action)							
EN 16615	high organic load		0,5 %		0,25 %		
activity against specific organisms: dermatophytes (skin fungi)							
EN 13624 (Trichophyton interdigitale)	high organic load			0,25 %			
additional activities and tests							
tuberculocidal (Mycobacter terrae)**				3 %	1 %		
Limited virucidal against enveloped viruses* (incl. HBV, HCV, HIV, human and animal influenza viruses)			0,5%				
Adenovirus*				0,5%	0,25%		
Vaccinia Virus*			0,5%				
BVDV*		0,25%					
Rotavirus*		0,25%					
Norovirus*						1%	
Polyomavirus SV 40*							0,5%
<b>VAH - Certificate</b>		<b>contact time</b>					
Surface disinfection in medical areas bactericidal and yeasticidal activity		5 min.	15 min.	30 min.	60 min.	120 min.	240 min.
	with mechanical action dirty conditions		0,5 %		0,25 %		

\*tested according DVV-method

\*\*DGHM-method

