

APESIN sanol

QUICK & EASY


Foaming sanitary disinfectant cleaner

- Descaling power
- Cleaning and disinfecting in one step
- Highly efficient

Product Profile

- Sanitary disinfectant cleaner with a broad activity against bacteria, yeast, viruses (limited spectrum) in medical, institutional, and public areas.
- APESIN sanol Quick & Easy dissolves quickly and easily mineral soiling like scale, lime soap deposits, and rust, as well as organic soiling, like grease and proteinic deposits.
- APESIN sanol Quick & Easy ensures effective disinfecting and cleaning results in one step.
- The high concentrated formulation minimizes transport weight, and leads to a lower inventory.
- The Quick & Easy portable dosing system offers accurate and safe dosage.
- By the elimination of ingredients, such as aldehyde, phenol, alcohol, and chlorine there is no odour nuisance.
- No aerosol thanks to innovation foam solution and therefore perfectly suitable in large-scale areas.

Area of application

- APESIN sanol Quick & Easy is suitable for all hard, and water resistant sanitary surfaces in medical, and institutional areas, e.g. public facilities, offices, gyms as well as swimming pools.
- The product can be used on POM/Polyoxymethylene Hostaform type 13031 XAS, type C 9021 GV 1/30.
- Also applicable on chrome tapware in good conditions, plus high-grade brass.
- Do not use on acid-sensitive materials, such as marble, copper, aluminium, nylon (polyamide), and PLEXIGLAS® / acrylic glass.

Application note

- Use separate trigger per APESIN product.
- Apply foam onto a cloth and completely wet the surface.
- Moisten the cloth sufficiently with APESIN sanol Quick & Easy in order to ensure appropriate wetting of the surface.
- Allow the solution to work for at least 5 minute.
- For detailed information see next page: Standards & spectrum of activity.

Ingredient declaration

Ingredients per 100 g: 9,9 g LACTIC ACID, < 5 % anionic surfactants, non-ionic surfactants, perfumes.

Sustainable development and production



Product safety, storage and environmental instructions

Safety: For professional use only. For further information, see Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

Storage: Store at room temperature in the original container. Keep from freezing.

Environmental: Only give completely emptied container to special waste disposer.



Application and Dosing



Use separate trigger for APESIN sanol Quick & Easy.



Dose in relation to application. Follow below instructions.



Empty water bottle after a working day.



Disinfection cleaning: Apply foam onto a cloth and completely moisten the surface. Allow solution to work according to exposure time stated on next page.



predefined dosing

Sales units:

Order No.: 715173

6 x 325 ml

Your Trusted Partner On Site:

pH-Value



APESIN sanol

**QUICK
& EASY**


Foaming sanitary disinfectant cleaner

Disinfection activity of APESIN sanol Quick & Easy

Surface Disinfection:		Exposure time				
Institutional, industrial & medical area		1 min.	5 min.	15 min.	30 min.	60 min.
Bactericidal						
EN 13727	High organic load		x			
EN 16615	High organic load		x			
Yeasticidal						
EN 13624	High organic load					x
EN 16615	High organic load		x			
		Exposure time				
		1 min.	5 min.	15 min.	30 min.	60 min.
Yeasticidal						
EN 14476 virucidal activity against enveloped viruses	High organic load		x			
EN 14476 (Vacciniavirus)	High organic load		x			

Quick & Easy system automatically ensures a correct dilution of min. 10 % (v/v) from the product concentrate.

*BPR: The Biocidal Products Regulation is valid from 2013. This regulation aims to improve the functioning of biocidal products market in the EU. For further information, please contact: info.professional@werner-mertz.com

