

LONGLIFE hospital

Disinfectant - resistant emulsion polish



Sustainability Information

Sustainable Contribution

- Developed to fulfill the Austrian Ecolabel
- Registered in the positive list of 'DIE UMWELTBERATUNG'
- Complete abandonment of toxic and harmful materials
- Outperforms other products at lower use concentration
- With surfactants of renewable plant-based origin



Quality & Performance

- More than 30 years of experience in the development of high-performance and integrally sustainable cleaning concepts.
- Our Green Care Professional products are produced exclusively in the EMAS-certified plants in Mainz and Hallein in compliance with the highest quality standards



Safe for humans & environment

- Our formulas and packaging are safe for humans and nature as well as biological and technical cycles when used properly
- LONGLIFE hospital is formulated from readily biodegradable raw materials



Circular cleaning products

- Production with intensive use of renewable raw materials and energy up to 100% from hydropower, solar and geothermal power
- To find out how much plastic, CO₂ and crude oil can be saved with Green Care Professional products, visit: get.wmprof.com



Sustainable profitability


- Our intelligent systems and efficient formulas enable precise and economical dosing as well as complete emptying of the packaging
- Comprehensive and consistently practiced philosophy of sustainable management



Comprehensive Service

- As a strategic partner, we offer comprehensive know-how transfer through holistic consulting and support as part of our Green Care Academy
- Customized hygiene plans and cleaning recommendations optimizing your cleaning procedures and thus your business

Product certificates



Company certificates of Management Systems



EMAS: Most demanding system for sustainable environmental



DIN ISO 14001: Internationally recognized environmental management system



DIN ISO 50001: Internationally recognized energy management system



DIN ISO 9001: Internationally recognized quality management system



EcoVadis: Provider of business sustainability ratings

ONGLIFE hospital

Disinfectant - resistant emulsion polish

Disinfectant-proof* | Long-lasting | Easy application



Product profile

- LONGLIFE hospital is a care dispersion that is resistant to disinfectants as well as foot marks, heel marks and other stresses.
- LONGLIFE hospital ensures optimum resistance to alcohol and disinfectants.
- The disinfectant-proof emulsion polish is time- and cost-saving and is suitable for areas with specific hygiene requirements.
- LONGLIFE hospital is recommended wherever an excellent look or durability of the protection film is desired.
- As a safety aspect LONGLIFE hospital is slip-resistant because of the contained powerful care components (according to DIN51131) even under strong frequentation.
- LONGLIFE hospital is highly resistant to heel marks and soils. It reduces resoiling and enables lasting cleanness.
- With its uniform look it is fast and simple in processing and shows an excellent flow further high-speed-polishing increases the durability.
- By the complete abandonment of hazardous ingredients such as zinc, TBEP, perfume, and fluorinated levelling agents, LONGLIFE hospital takes on the responsibility for future generations.

Area of application

- LONGLIFE hospital is very suitable for the maintenance of water-resistant floorings especially in health facilities, for example hospitals, doctor's surgeries, retirement and nursing homes, rehab centers etc.
- When properly applied, LONGLIFE hospital creates a hard-wearing protective film with optimum resistance to disinfectants. In areas of dosing dispensers, dripping of disinfectants onto the floor can never be completely avoided. Due to the good alcohol resistance, the care film remains largely unchanged even after the drying of alcoholic products.
- Do not apply on unsealed wood. Before using the product on elastomer coverings please install a testing area.

Ingredients, handling, storage and environmental instructions

For the ingredients list, please see the Safety Data Sheet. Product strictly reserve for a professional use. For all information regarding handling, storage conditions and disposal, please see the Safety Data Sheet.

Product variants

716038 - 2x 5 l


pH value -



Application and Dosing




Dose in relation to application and degree of soil. Follow below instructions.



concentrated


Floor preparation:

Strip the floor, remove all dirt and old dispersion layers.




Initial treatment:

Apply 2-3 thin and uniform layers. Drying time between the applications at least 30min.



Large areas:


Work in sections and apply wet in wet. Usage: approx. 2l / 100 m².



200 ml/l

Coating refreshment:

Remove slight scratches using the spray-cleaner-method and a highspeed machine with red pad.



After awaiting the fair drying time (one day) polish with Ultra High-Speed-Method (>1000 U).

Learn more at: www.green-care-academy.com