

APESIN handactive F

Liquid alcoholic instant hand disinfectant

Hygienic, surgical handrub | Ready-to-use



Technical Information



Product profile

- APESIN handactive F is a ready-to-use perfume- and dye-free hand disinfectant that ensures fast-acting hygienic and surgical hand disinfection.
- It contains moisturizing and skin-caring additives.
- APESIN handactive F is effective against bacteria, mycobacteria, yeasts and enveloped viruses.

Area of application

- **Hygienic handrub:** Indoor - hospitals and other health care institutions, ambulances, surgeries, nursing homes (including home-care of patients). Hospital canteens, large kitchens, pharmaceutical industries, production sites, laboratories: hygienic handrub onto visibly clean and dry hands.
- **Surgical handrub:** Hospitals and other health care institutions: surgical handrub onto visibly clean and dry hands and forearms.
- Authorization number: EU-0027673-000

Application note

• PT 1 - Human hygiene

Number and timing of application: There is no restriction in the number and timing of applications. No safety intervals need to be considered between the application phases. The product may be used at any time and as often as required.

- Hygienic handrub:

Application Rate: Dosage: At least 3 ml. Contact time: 30s.

- Hygienic handrub onto visibly clean and dry hands.

- Surgical handrub:

Application Rate: Dosage: Rub sufficient amount in portions of 3 ml. Contact time: 90 s.

Sales units:

1114909 - 10 x 1 l

1116199 - 10 x 500 ml



Application and Dosing



Hygienic handrub

Apply at least 3 ml of product onto clean and dry hands and ensure complete wetting of the hands during the whole exposure time. For more information, see exposure times stated on next page.



Surgical handrub

For surgical handrub use as many portions of 3 ml as necessary to keep hands wet for 90 seconds.

APESIN handactive F

Hand disinfectant		
hygienic handrub		
target organisms	Dosage / Dilution	contact time
Bacteria	At least 3 ml/ Ready-to-use product	30 s
Mycobacteria		
Yeasts		
Enveloped viruses		
surgical handrub		
target organisms	Dosage / Dilution	contact time
Bacteria	Sufficient amount in portions of 3 ml/Ready-to-use product	90 s
Mycobacteria		
Yeasts		
Enveloped viruses		

The results have been substantiated to the competent authority by testing the following standards: EN 13727, EN 13624, EN 1500, EN 12791, EN 14348, EN 14476

VAH -Certificate	Dosage / Dilution	contact time
hygienic handrub	2 x 3 ml Ready-to-use product	30 s
surgical handrub	Sufficient amount in portions of 3 ml/Ready-to-use product	90 s

Ingredient declaration

Ingredients per 100g: 45,0 g propan-2-ol, 30,0 g propan-1-ol.

Handling, storage and environmental instructions

Product strictly reserved for a professional use. For all information regarding handling, storage conditions and disposal, please see the Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

Company certificates of Management Systems

