

APESIN daily



Liquid disinfecting cleaner

- Versatile use
- Surface disinfection for industrial and institutional areas

Product profile

- Cleaning and disinfection in one step of floors and surfaces in industrial and institutional areas
- Broad biocidal spectrum with biocidal, yeasticidal and virucidal activity
- Virucidal activity against enveloped viruses such as e.g. Vaccinia Virus, Influenza viruses, Corona viruses (incl. SARS-CoV & MERS-CoV), Hepatitis, Ebola Virus, etc.

Application

- For the correct dosage, please follow the below instructions.
- Apply diluted cleaning and disinfection solution with a mop or sponge and let solution act for the required disinfection time.
- Apply after every working shift or when required during the working processes.
- For the usage in food processing areas with surfaces in direct contact with food rinsing with water is required prior to the re-use of the surfaces.

Ingredient declaration

Ingredients per 100g: 4,5g DIDECYLDIMONIUM CHLORIDE, 5 - <15% Non-ionic surfactants

Sustainable development and production



Product safety, storage and environmental instructions

Safety: For professional use only. For further information, see Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

Storage: Store at room temperature in the original container.

Environmental: Only give completely emptied packaging to special waste disposer.

Application & dosing



Dose in relation to application and degree of soil. Follow below instructions.



Floor cleaning: Wipe floor wet or damp.



Ideal for application on surfaces in kitchen areas.



Applicable in sanitary areas.



Floor

Sales units:

Bestell-Nr.: 713957 2 x 5 L

Your trusted partner on site:

pH-value:



APESIN daily

Liquid disinfecting cleaner

Surface disinfection:		Exposure time				
Food area, industrial & institutional area		1 min.	5 min.	15 min.	30 min.	60 min.
Bactericidal activity						
EN 1276	high organic load		3%			
EN 13697	high organic load		3%			
Yeasticidal activity						
EN 1650	high organic load			2,5%		
EN 13697	high organic load			2%		
Virucidal activity						
EN 14476	low organic load		1,0%			
virucidal activity against enveloped viruses						
EN 14476 individual test organisms						
Vaccinia virus	low organic load		1,0%			



APESIN daily

Liquid disinfecting cleaner

Surface disinfection:		Exposure time				
Food area, industrial & institutional area		1 min.	5 min.	15 min.	30 min.	60 min.
Bactericidal + yeasticidal activity						
EN 1276 + EN 1650 + EN 13697	high organic load			2,5%		
Virucidal activity						
EN 14476 effective against enveloped viruses (incl. HIV, HBV, HCV)	low organic load		1,0%			
EN 14476 individual test organisms						
EN 14476 Vaccinia Virus	low organic load		1,0%			

